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90mg/3mL Sodium Hyaluronate

## DESCRIPTION OF DEVICE

Hyaluronic acid is distributed throughout all body tissues and is found in high concentrations in tissues such as the vitreous, synovial fluid, umbilical cord and dermis. Hyaluronic acid acts as a tissue lubricant and is thought to play an important role in regulating interactions between adjacent tissues. It also acts as a viscoelastic support that keeps tissues separate. Different hyaluronic acid preparations may have different molecular weights but are thought to have the same chemical structure. The product does not interfere with normal wound healing processes. The device is invasive and implant device.

## PURPOSE OF USE and INDICATIONS

HyaRelief™ 3.0% is intended for the treatment of symptoms of human joint disorders such as osteoarthritis.

HyaRelief™ 3.0% is indicated for viscoelastic replenishment or replacement of synovial fluid in human joints.

## CONTRAINDICATIONS

Known sensitivity to any of the ingredients contained in product,  
Pre-existing infections of the skin in the region of the intended injection site,  
Known infection of the knee joint.  
Known systemic bleeding disorders.  
Do not use on infected and inflamed joints.

HyaRelief™ 3.0% may contain trace amounts of gram positive bacterial proteins and are contraindicated for patients with a history of such allergy.

## ADVERSE REACTIONS AND COMPLICATIONS

Infection, arthralgia (knee pain), arthrosis, joint (knee) disorder, joint (knee) swelling, joint (knee) effusion, joint (knee) stiffness, pain in limb, tendonitis, paraesthesia, phlebitis, pruritus, injection site erythema, injection site edema, injection site pain, injection site reaction, arthropathy, baker's cyst, bursitis, localized osteoarthritis, aggravated osteoarthritis, immune response.

## PATIENT PROFILE

HyaRelief™ 3.0% is recommended for patients with Osteoarthritis and Joint Diseases. Product is not used in pregnant or breast-feeding woman, and children. The patients who under chemotherapy shall be considered by orthopaedist.

## USER PROFILE

Intra-articular Injection should be applied by Orthopedists.

## STORAGE CONDITIONS

The products are stored at 2-25 °C and low humidity (40-60% RH) before being shipped to the customer. Excessive storage conditions such as exposure to direct sunlight and/or storage on or near radiators may destroy the product.  
Do not freeze the product.

## WARNINGS AND PRECAUTIONS


### PRODUCTS INTENDED FOR SINGLE USE MUST NOT BE RE-USED.

- Re-use or reprocessing (e.g., cleaning and re-sterilization) may compromise the structural integrity of the device and/or lead to device failure which may result in patient injury, illness or death.
- Furthermore, reuse or reprocessing of single-use devices may create a risk of contamination e.g. due to the transmission of infectious material from one patient to another. This could result in injury or death of the patient or user.
- Contaminated implants must not be reprocessed. Any product that has been contaminated by blood, tissue, and/or bodily fluids/matter should never be used again and should be disposed of according to hospital protocol.
- Only medical professionals trained in accepted injection techniques for delivering agents to intra-articular synovial joint spaces should inject sodium hyaluronate for this application.
- Injection should be administered according to the severity of the disease, the volume of the application area and the decision of the physician (orthopedist).
- If pain increases during the injection procedure, the injection should be stopped, and the needle withdrawn.
- Patients experiencing abnormal sequelae (tissue defect from a disease) to the administration of product should consult with the physician immediately.


- The device should be used with 18-22G x½ needle.
- Syringes, needles, and any remaining product should be disposed of in accordance with the legal requirements for medical waste status in the respective countries.
- The physician must follow the rules of the technique of intra-articular injection with suitable angle and site of the joint space. Do not inject into glandular tissue.
- Do not concomitantly use disinfectants containing quaternary ammonium salts for skin preparation because sodium hyaluronate can precipitate in their presence.
- Inject into the synovial space only. Do not inject by intravascular route.

## HOW TO USE

- The product should be administered according to the frequency determined by the physician, according to the recommended frequency decision, such as 6-8-12 months.
- In case of joint effusion, aspiration, rest, ice pack application and/or intra-articular cortisone administration is recommended to reduce effusion. After 2-3 days, product application is started.
- The contents of product and the outer surface of the pre-filled syringe are sterile until the package is opened.
- Take the syringe out of the package.
- Remove air bubble (if any) before injection.
- Remove the Luer-lock hood.

- Attach a suitable cannula (18-22G x½) carefully by rotating it.

## PACKAGE CONTENT

The package is included 1 pre-filled syringe.



## STERILIZATION and SHELF-LIFE

The device is sterilized by autoclave sterilization.

The shelf-life is 2 years.

Release Date: 12.02.2021 Revision No: 04 Revision Date: 01.07.2024

## SYMBOLS

	CAUTION		DO NOT USE IF PACKAGE IS DAMAGED
	SEE INSTRUCTIONS FOR USE		EXPIRATION DATE
	DO NOT RE-STERILE		KEEP AWAY FROM SUNLIGHT
	LOT NUMBER		TEMPERATURE LIMITATION 2°C - 25°C
	PRODUCTION DATE		STERILIZED USING MOIST HEAT
	CONFORMITY EUROPEAN AND NOTIFIED BODY NUMBER		STERILE FLUID PATH
	MEDICAL DEVICE		HUMIDITY LIMITATION %40 - %60
	DRY KEEP		NON-PYROGENIC
	DO NOT RE-USE		REFERENCE NUMBER
	MANUFACTURER		UNIQUE DEVICE IDENTIFIER